

Biochip Course

Optical Detection for Biosensing

Marco Carminati

a.y. 2023/2024

April 8th 2024




Detection Approaches

Extremely **thin** sensitive layers, in liquid environment



Detection techniques sensitive to planar surface changes



Everything happens at the **interface!**





Optical microscopy is the preferred detection tool in Biology
non invasive, relatively high power, molecularly selective

1. Fluorescence: GFP, immunostaining
2. Scanning probe microscopy
3. Sources, filters, detectors
4. Solid-state integrated detection: Silicon Photonics



Fluorescence

Light emission by optically excited molecules.




A large energetic separation (*Stokes shift*) makes easier the detection

Bleaching limits the number of excitation cycles



Fluorescence Microscope



Optical microscopy is the workhorse of laboratory-scale Biology research.


Since cells are transparent, fluorescence microscopy is used to detect sub-cellular entities.

It can report the presence of a given molecule bound to the fluorescent molecule.





Optical Marker: Green Fluorescent Protein

Fluorescent protein of natural origin
used as a photo-activable
label of genetic expression



Aequorea Victoria



Different Molecules and Filters

Engineered more molecules to provide more colors.

The microscope has a set of corresponding filters (R,G,B)

Excitation Illumination Balancer Filter Configuration and Spectral Profiles






Figure 1






Examples



Bright Field



Fluorescence



Combined





Selective staining: penetrate the cell membrane and bind to specific parts:

- Blue: nucleus
- Green: cytoskeleton
- Red: mitochondria



Molecular Markers Detection: Immunostaining


The optical label is attached to an antibody that **specifically** binds with a target molecule (antigen), whose presence can be optically detected.





Enzyme-Linked Immunosorbent Assay


Standard biochemical test: the color change is proportional to the quantity of the enzyme






Types of Assays

Sandwich



Optical (fluorescent or visual) signal is proportional to the quantity of bound analytes

Competitive




Signal decreases proportionally to the analytes replacing the bound markers




The Limits of the Optical Microscope



The spatial resolution is limited by diffraction:



Single
fluorophore
(4 nm) is a
point source





$$r = \frac{0.5 \lambda}{NA} = \frac{0.5 \lambda}{n \sin(\theta)}$$




Beyond The Limits of Diffraction

Surfaces can be “imaged” below the optical resolution limit, by scanning them with a sharp microfabricated tip, interacting with the surface atoms via short-range **forces** whose intensity is **spatially mapped**




Scanning Probe Microscopy







A Family of Techniques



Binnig and Rohrer,


Nobel Prize in Physics, 1986

The Atomic Force Microscope



... like reading Braille!




Optical Detection in Biochips





Light Sources

- **Lamps**
 - High power
 - High dissipation
 - Wide spectrum (filter)
- **LEDs**
 - Solid-state
 - High efficiency
 - Limited power (growing)
- **Lasers**
 - Monochromatic
 - Solid-state
 - cost/size, few λ








Absorption Optical Filters

Colored (pigmented) films that absorb (i.e. convert into heat) some wavelengths

- Large bandwidth
- Shallow slopes



Broad Band Absorption Filter Spectra





Interference Filters

Stack of thin layers with alternate refractive index

- Controllable properties (nm deposition)
- Sharp slopes

Reflection and Transmission by Interference Filters




Figure 6

Refraction and Reflection in Optical Coatings





Figure 4



Effect of Cavity Number on Pass





Light Detectors

- Photoresistors
- **Photodiodes:**
 - PN
 - PIN
- Avalanche photodiodes:
 - Photomultipliers (PMT)
 - APDs
 - SPADs
 - SiPMs
- Arrays (image sensors):
 - CCD
 - CMOS





Detection Approaches for Biosensing

With external stimulation:

- Light absorption (light extinction)
- Fluorescence (energy transfer)
- Evanescent field interactions (light propagation)
- *Surface Plasmon Resonance*: sophisticated technique very sensitive to the metal surface properties.


No external source excitation:

- *Chemiluminescence*: some biochemical reactions produce light as a product




Spectroscopy

Each compound has an absorbance spectrum (fingerprint):




Spectrophotometer:



Requires a dispersive medium (prism, grating) and a broadband detector (or an array)



At a single wavelength:



Lambert-Beer Law

Absorption:

$$A = \log_{10} \left(\frac{I_0}{I} \right) = \log_{10} (\%T^{-1}) = \epsilon l c$$

$\alpha = A/l$ = absorbance

l = length


c = concentration of molecules

ϵ = molar absorbance



Example: Pulse Oximetry

$$\begin{aligned}\alpha(\lambda_1) &= \varepsilon_{\text{ox}}(\lambda_1)[M_{\text{ox}}] + \varepsilon_{\text{deox}}(\lambda_1)[M_{\text{deox}}] + \varepsilon_{\text{meth}}(\lambda_1)[M_{\text{meth}}] \\ \alpha(\lambda_2) &= \varepsilon_{\text{ox}}(\lambda_2)[M_{\text{ox}}] + \varepsilon_{\text{deox}}(\lambda_2)[M_{\text{deox}}] + \varepsilon_{\text{meth}}(\lambda_2)[M_{\text{meth}}] \\ \alpha(\lambda_3) &= \varepsilon_{\text{ox}}(\lambda_3)[M_{\text{ox}}] + \varepsilon_{\text{deox}}(\lambda_3)[M_{\text{deox}}] + \varepsilon_{\text{meth}}(\lambda_3)[M_{\text{meth}}]\end{aligned}$$




External monitor of O_2 saturation in blood




How to Handle Light on Chip: Silicon Photonics


Monolithic integration of sources, detectors, waveguides, switches, filters etc... on a single solid-state active (or passive) platform (Si, InP, Ge...)

Applications: telecom (long and short distance), processing and also **biosensing**





Confined propagating EM field




1. Total internal reflection condition

$$n_1 \sin \theta_1 = n_2 \sin \theta_2$$

$$\theta_1 > \arcsin\left(\frac{n_2}{n_1}\right)$$

2. Phase coherence wavefronts




$$\frac{2\pi}{\lambda} n_w d \cos \theta + 2\phi = N\pi$$


Integrated Waveguide Types

can be air


Buried





Ridge



Rib



Diffused






Materials:

- Core: Si, Ge, InP
- Cladding: SiO_2 , air...



Evanescent Field




At total internal reflection,

Exponential decay of field intensity

Propagating light is **affected by adsorbed molecules**





Waveguide as Biosensor





The Mach-Zehnder Interferometer


Useful to convert phase difference into intensity difference:



Light splits and combine after two different paths, whose optical length can be controlled (modulated)




Interferometer Branch Used as Biosensor







Implementation





Optical Ring Resonators

Take advantage of resonance



Resonance condition: $n_w L = N\lambda$

$$FSR = \frac{c}{n_w L}$$

Q depends on loss and coupling





Free Spectral Range

Biosensors Ring Resonators

The presence of the analyte changes n






Change in the resonance frequency





A Transparent In-Line Light Power Monitor



Invented at Polimi!

$$\lambda[\mu m] = \frac{1.24}{E_{gap}[eV]}$$



$$C_A = \epsilon_{SiO_2} \frac{L_e w}{t_{ox}} \quad G = \sigma_{Si} \frac{h w}{L_e}$$

Although Si is transparent at 1550nm, the defects at the surface produce inter-gap energy states that allow photon absorption

Insulated far top metal electrodes do not perturb the optical field



A Transparent In-Line Light Power Monitor



Local power intensity can be monitored by a contactless AC measurement of the conductivity of the Si waveguide

$$f_e > \frac{1}{2\pi C_A / 2R_{WG}} \quad R_{WG} = \frac{1}{G}$$




Key Take-Home Messages

- Fluorescent microscopy leverages fluorescent molecular markers (first GFP) that can be genetically manipulated.
- In a biochip the optical system comprises: a source (LED, laser, lamp) and a photodetector.
- To map **surfaces** below the optical diffraction limit, scanning probe techniques (AFM in particular) can be used.
- Integrated photonic waveguides (arranged in interferometer or resonators devices) whose cladding is functionalized with receptors can be used to detect molecular binding. The quantitative relation between the analyte concentration and Δn is critical and need to be elucidated for each sensor.



References



Silvano Donati, Photodetectors: Devices, Circuits, and Applications, Prentice Hall.

IEEE JOURNAL OF SELECTED TOPICS IN QUANTUM ELECTRONICS, VOL. 20, NO. 4, JULY/AUGUST 2014

8201710

Non-Invasive On-Chip Light Observation by Contactless Waveguide Conductivity Monitoring

Francesco Morichetti, Stefano Grillanda, Marco Carminati, *Member, IEEE*, Giorgio Ferrari, *Member, IEEE*, Marco Sampietro, Michael J. Strain, Marc Sorel, and Andrea Melloni